



Lebrush Lab

Lebrush Lab

Stability and Innovation—Leading the Way to New Possibilities.



About Us

LeBrush Lab was founded with a clear mission: to address the challenges posed by the high and often fluctuating cost of Kolinsky hair, as well as its inconsistent supply. As the market began seeking alternatives, we recognized a unique opportunity for innovation. After years of research and development, we discovered a natural material that provides a stable and economical alternative to Kolinsky hair.

Our journey has been built not only on advances in material science but also through close collaboration with dental technicians and professionals. Through continuous testing and refinement, we ensure that every brush meets the rigorous standards required in practice. Feedback from experts has guided us in optimizing our product, achieving ideal levels of flexibility, durability, and smoothness to fully meet the precise needs of dental technicians in porcelain work.





Harmony Porcelain Brush

Perfect Harmony, Enduring Quality.



Harmony Porcelain Brush



Harmony Porcelain Brush

LLVS04-1

\varnothing 3.0mm
 \leftrightarrow 14.0mm

LLVS06-1

\varnothing 4.4mm
 \leftrightarrow 17.0mm

LLVL06-1

\varnothing 4.4mm
 \leftrightarrow 17.0mm

LLVS08-1

\varnothing 5.5mm
 \leftrightarrow 20.0mm

LLVL08-1

\varnothing 5.5mm
 \leftrightarrow 20.0mm



Harmony Porcelain Brush



Harmony Porcelain Brush

LLVS04

Ø 3.0mm
↪ 14.0mm

LLVS06

Ø 4.4mm
↪ 17.0mm

LLVL06

Ø 4.4mm
↪ 17.0mm

LLVS08

Ø 5.5mm
↪ 20.0mm

LLVL08


Ø 5.5mm
↪ 20.0mm

Harmony Porcelain Brush




L08

The thick brush tip is suitable for applying porcelain on large surfaces, improving efficiency while enhancing the strength and aesthetics of the dental prosthesis.




S08

Ensures a natural transition and perfect fusion of porcelain layers, with precise powder layering in hard-to-reach areas, restoring the tooth's anatomy realistically.



L06

Balances precision and coverage, ideal for crown surfaces, providing even porcelain thickness with a natural color transition.



S06

The long-tipped design aids in shaping complex tooth structures, ensuring precise porcelain layering and a smooth, natural form.





Harmony Porcelain Brush



LLVS000T



LLVS01T



LLVS02T



LLVS04T



Fibertec
Wood Handle

LLVS000T

Ø 0.5mm
↪ 6.0mm

LLVS01T

Ø 0.8mm
↪ 6.8mm

LLVS02T

Ø 1.2mm
↪ 9.0mm

LLVS04T

Ø 1.5mm
↪ 10.0mm

LLVS000T

Ø 3.0mm
↪ 14.0mm



Harmony Porcelain Brush



LLVS06T



LLVL06T



LLVS08T



LLVL08T

Brush Tip
(2 pcs pack)

LLVS06T

Ø 4.4mm
↪ 17.0mm

LLVS08T

Ø 5.5mm
↪ 20.0mm

LLVL06T

Ø 4.4mm
↪ 17.0mm

LLVL08T

Ø 5.5mm
↪ 20.0mm

Fibertec Porcelain Brush

Engineered Fibers, Built to Last.





Fibertec Porcelain Brush



Fibertec Porcelain Brush

LLSS02

\varnothing 1.5mm
 \downarrow 10.0mm

LLSS04

\varnothing 3.0mm
 \downarrow 14.0mm

LLSS06

\varnothing 4.4mm
 \downarrow 17.0mm

LLSL06

\varnothing 4.4mm
 \downarrow 17.0mm

LLSS08

\varnothing 5.5mm
 \downarrow 20.0mm

LLSL08

\varnothing 5.5mm
 \downarrow 20.0mm



Fibertec Porcelain Brush



LLSSF01



LLSS000



LLSS0



LLSS01

Fibertec Porcelain Brush

LLSSF01

∅ 3.7mm
↪ 7.0mm

LLSS000

∅ 0.8mm
↪ 6.0mm

LLSS0


∅ 0.8mm
↪ 8.0mm

LLSS01

∅ 1.2mm
↪ 9.0mm



Fibertce - Wood Handle



Fibertec Wood Handle

LLNS02

\varnothing 1.5mm
 \mapsto 10.0mm

LLNS04

\varnothing 3.0mm
 \mapsto 14.0mm

LLNS06

\varnothing 4.4mm
 \mapsto 17.0mm

LLNL06

\varnothing 4.4mm
 \mapsto 17.0mm

LLNS08

\varnothing 5.5mm
 \mapsto 20.0mm

LLNL08

\varnothing 5.5mm
 \mapsto 20.0mm



Fibertec - Wood Handle



LLNSF01



LLNS000



LLNS0



LLNS01

Fibertec Wood Handle

LLNSF01

∅ 3.7mm
↪ 7.0mm

LLNS000

∅ 0.5mm
↪ 6.0mm

LLNS0


∅ 0.8mm
↪ 8.0mm

LLNS01

∅ 1.2mm
↪ 9.0mm



Harmony - Wood Handle



Harmony Wood Handle

LLTS02

\varnothing 1.5mm
 \downarrow 10.0mm

LLTS04

\varnothing 3.0mm
 \downarrow 14.0mm

LLTS06

\varnothing 4.4mm
 \downarrow 17.0mm

LLTL06

\varnothing 4.4mm
 \downarrow 17.0mm

LLTS08

\varnothing 5.5mm
 \downarrow 20.0mm

LLTL08

\varnothing 5.5mm
 \downarrow 20.0mm



Fibertec - Wood Handle



LLTSF01



LLTS000



LLTS0



LLTS01

Fibertec Wood Handle

LLTSF01

∅ 3.7mm
↪ 7.0mm

LLTS000

∅ 0.5mm
↪ 6.0mm

LLTS0


∅ 0.8mm
↪ 8.0mm

LLTS01

∅ 1.2mm
↪ 9.0mm



SET



Set Brush: #8 complete brush,
with brush tips size: #000 , #01, #04

SET

LLVS-1 (Red)

LLVS-2 (Blue)




SET

LLSS-1



LLNS-1



SET

LLSS-1 (Acrylic Handlen)

LLNS-1 (Wood Handle)


#Set Brush : Complete brush#8, #000,
#01, #04



Lebrush Lab

Lebrush Lab

Stability and Innovation—Leading the Way to New Possibilities.



Instagram



WhatsApp



WeChat